


PORTER...
the name to *trust* in Nitrous

MXR-D DIGITAL FLOWMETER



PORTER MXR-D Precision Nitrous Oxide Sedation
with visual confirmation of real flow

PORTER MXR-D

The first **electronic** flowmeter providing the assurance of **visual flow confirmation**, as well as **precise digital readouts** of the actual flow and gas mixtures.



Fits or retrofits easily into existing or new cut-out templates



Needle valve technology

- Instantaneous flow of gases at the turn of a knob.
- % N₂O adjustment knob: Controls the % concentration of N₂O. (70% maximum N₂O)

Digital display readout.

- Normal flow display shows the % N₂O on the left side and total flow in L/min on the right side. Additional user conditions such as "O₂ flush" and "Low O₂ flow" will also display.

On/Off switch

- Provides a manual on/off switch with a dual function as the electrical on/off for the digital display.

Total flow adjustment knob

- Controls the combined flow of O₂ and N₂O (or O₂ flow only when the % N₂O knob is set to zero).

Flushmount design contributes to:


- ✓ Ease of use
- ✓ Frequency of use
- ✓ Patient satisfaction
- ✓ Higher case acceptance
- ✓ More "billable" procedures, and
- ✓ Increased sales and gross profits

Oxygen flush button

- Easily accessible button, when pressed, provides high flow of oxygen.

A leap forward in precision and safety from Porter, the leader in N₂O sedation.

- ✓ Porter's proprietary needle valves deliver instantaneous flow vs. the delay in the electronic valving of digital flowmeters.
- ✓ Because of Porter's control valve design, electrical surges or power outages have no effect on the flow of gas.
- ✓ Audible alarms warn of any interrupted flow conditions. A minimum 30% O₂ flow is always maintained.
- ✓ The MXR-D mounts easily in all cabinetry.
- ✓ Safe, low voltage design (12v DC).
- ✓



Digital Flowmeter Comparison Chart

	Manufacturer		
	Porter MXR-D	Digital Flowmeter "A"	Digital Flowmeter "M"
Digital Display	YES	YES	YES
Calibrated Flow Tubes	YES	NO	NO
Flow control: 0-6 Lpm (time)	< 1 sec	> 6 sec	> 6 sec
Power interruption flow protection	YES	NO	NO
Real time flow assurance	YES	NO	NO
O ₂ quick-connect on meter	YES	NO	NO
Power to the meter	12v DC	110v AC	110v AC

O₂ quick connect

- Facilitates connection of resuscitation equipment to the central system for emergency oxygen. Comes standard with Porter, an upgrade with other meters.



Calibrated Real Flow

Porter N₂O and O₂ indicator tubes

- Calibrated flow tubes provide continuous flow in case of power interruptions.



Unsursed Scavenger Efficiency



Patented AVS: Automatic Vacuum Switch

"The Porter Cabinet Mounted Flowmeters are so convenient to use, I find myself using nitrous more with my patients. While I can't quantify the percentage nitrous has contributed to our office's 'full appointment book', I do know a cooperative patient allows me to do my best dentistry. Our patients have a relaxing experience instead of a fearful encounter."


Jon Scharf, DDS - Exton, PA

The Porter MXR-D Digital Flowmeter, an integral part of the Porter Nitrous Oxide Sedation System.

- Better for your patients.**
- Better for your practice.**
- Better for your bottom line.**

The Porter Nitrous Oxide Sedation System

Better for your patient. Better for your practice. Better for your bottom line



PORTER

Parker Hannifin Corporation
Porter Instrument Division
 245 Township Line Road
 P.O. Box 907
 Hatfield, PA 19440-0907 USA
 Phone (215) 723-4000
 Fax (215) 723-2199
www.porterinstrument.com